

A close-up photograph of a woman with long dark hair looking down at a baby with blonde hair. The baby is wearing a pink knitted sweater. The word 'ora' is overlaid in a large, white, serif font.

ora

Non-Invasive Endometrial Receptivity Test

ora[™] A Revolution in Endometrial Receptivity Testing

What is *ora*[™]

ora[™] is the world's first non-invasive endometrial receptivity test for identifying a patient's optimal window of implantation (WOI). It assesses a combination of microRNA (miRNA) biomarkers in the bloodstream and physiological conditions to determine the status of a patient's endometrium, providing information that can be used to optimize the timing of implantation.

Who is *ora*[™] for


If any of the following situations apply to your patients, *ora*[™] may be able to help:

- A history of implantation failure or miscarriage
- Few remaining high-quality embryos
- A lower or higher BMI
- Age 35 or older*

* While all ages have the possibility of having a displaced WOI, recent data has shown that this is an increased possibility for patients 35 or older.

Window of Implantation

The average WOI is between days 19 and 21 of the menstrual cycle. However, this can vary among individuals. Among women who have trouble conceiving, around 30% have been found to have a displaced WOI – one that occurs earlier or later than average.



Why choose *ora*™ ?

ora™ uses novel miRNA biomarkers present in the blood to accurately identify a patient's WOI, all with simple blood samples taken on day 4 and day 5 of a mock cycle*. This unique endometrial receptivity test delivers fast and reliable results, helping increase the successful pregnancy rate for IVF patients with a history of implantation failure.

* We recommend drawing blood on both day 4 and day 5 to account for the possibility of a shortened WOI post-receptive result. However healthcare providers may choose to do only a single blood draw on day 5.



NON-INVASIVE

ora™ analyzes miRNA biomarkers in the blood, removing the need for an invasive and uncomfortable endometrial biopsy procedure.



SINGLE MOCK CYCLE

By drawing blood on day 4 and day 5, *ora*™ is able to account for both average and shortened post-receptive windows, meaning patients only require one mock cycle regardless of what stage their endometrium is in.



HIGH ACCURACY


ora™ has been shown to have > 95% accuracy in predicting endometrial receptivity, offering a stable solution with a more comfortable testing procedure. Further, the need to re-test on account of inconclusive results or invalid/insufficient RNA occurs in <1% of cases.



COMPREHENSIVE ANALYSIS

ora™ analyzes close to 300 miRNA biomarkers that target over 1,000 endometrial receptivity-related genes to accurately identify the optimal time for embryo transfer.

Why use miRNA biomarkers?




MiRNAs have been shown to have over 90% accuracy in predicting displaced window of implantation as cause for implantation failure (1).

- ✓ **microRNAs (miRNAs) are enclosed and protected by proteins, making them more stable than mRNAs, facilitating analysis of lower-quality endometrial tissue samples.**
- ✓ **miRNA exhibits a high correlation with the protein level, as miRNAs regulate mRNAs to suppress protein translation and/or induce mRNA degradation.**
- ✓ **Numerous scientific publications have indicated that miRNAs play a role in the embryo implantation process, including regulating the endometrium to prepare for an embryo.**

miRNAs in tissue and blood

Endometrial receptivity tests were originally performed using endometrial tissue samples. microRNA (miRNA) biomarkers in the bloodstream can identify those same endometrial conditions.



These miRNAs in the bloodstream have been found to regulate a number of immune mechanisms during pregnancy, and are critical physiological factors for cell growth and angiogenesis. Through the regulation of these mechanisms, miRNAs can affect the endometrium's environment, growth process, and, by extension, pregnancy results.

Inti Labs was able to identify blood-based miRNA biomarkers that accurately reflect the endometrium's status, similar to the traditional method of analyzing endometrial tissue samples, but without the invasive biopsy needed for sample collection.

What can *ora*[™] tell you?

ora[™] uses these miRNA biomarkers to accurately analyze the endometrium's status, providing healthcare providers with a more precise recommendation for the optimal time for embryo transfer.



Receptive

The day 5 blood draw was found to be within the window of implantation (WOI), and is the optimal time for embryo transfer.



Pre-Receptive

The blood was drawn before the WOI, and the endometrium was not yet ready for embryo transfer. It is recommended to delay the embryo transfer cycle by 24 hours.



Post-Receptive (short)

The blood was drawn after the WOI, and the endometrium had passed the optimal time for embryo implantation. Based on the day 4 and day 5 blood draw results, it is recommended to move the embryo transfer time forward by **12 hours**.



Post-Receptive (average)

The blood was drawn after the WOI, and the endometrium has passed the optimal time for embryo implantation. Based on the day 4 and day 5 blood draw results, it is recommended to move the embryo transfer time forward by **24 hours**.

Inconclusive:

The resulting data could not be analyzed by *ora*[™]'s algorithm. This may be due to an exceptionally low-quality sample or complicated physiological conditions that are affecting the sample. An *ora*[™] representative will follow up with the healthcare provider to discuss performing the blood draw(s) again.

Invalid/insufficient RNA:

Results cannot be obtained due to the low quality or low concentration of the blood sample. The blood draw should be performed again to obtain a higher quality or higher concentration sample. An *ora*[™] representative will follow up with the healthcare provider to discuss re-performing the blood draw.

References

For Tissue miRNA

Chen CH, Lu F, Yang WJ, et al. A novel platform for discovery of differentially expressed microRNAs in patients with repeated implantation failure. *Fertil Steril*. 2021;116(1):181-188. doi:10.1016/j.fertnstert.2021.01.055

Chen CH, Lu F, Yang WJ, et al. Development of a Novel Endometrial Signature Based on Endometrial microRNA for Determining the Optimal Timing for Embryo Transfer. *Biomedicine*. 2024;12(3):700. Published 2024 Mar 21. doi:10.3390/biomedicine12030700

For Blood miRNA

Chen MJ, Hsu A, Lin PY, et al. Development of a Predictive Model for Optimization of Embryo Transfer Timing Using Blood-Based microRNA Expression Profile. *Int J Mol Sci*. 2023;25(1):76. Published 2023 Dec 20. doi:10.3390/ijms25010076



06 552 552 1